

MATERIAL SAFETY DATA SHEET

SECTION I

MANUFACTURER:

Co-Oral-ite Dental Mfg Co, Inc.
3340 Sunrise Blvd
Rancho Cordova, CA 95742

TELEPHONE:

(916) 638-8282

PRODUCT:

Teets Denture Material Powder - Cold Cure
~~WHITE DENTURE MATERIAL POWDER - COLD CURE (Plasticized~~
Methacrylate)

FORMULA: $-(CH_2C(CH_3)COOCH_3)-X$ + Dialkyl Phthalate

SECTION II: PRINCIPAL HAZARDOUS INGREDIENTS

	%	TLV
Dialkyl Phthalate	<15	5 ppm

SECTION III: PHYSICAL DATA

EVAPORATION RATE: Nil

VAPOR DENSITY: Not Applicable

SOLUBILITY IN WATER: Insoluble

PERCENT VOLATILE: <5%

APPROXIMATE BOILING RANGE: Not Applicable

DENSITY: Not Applicable

APPEARANCE AND ODOR: Fine powder. Faint odor in bulk.

SECTION IV: FIRE & EXPLOSION DATA

FLASH POINT (METHOD): Flash ignition approx. 304C (TCC)

APPROX. FLAMMABLE LIMITS: Not applicable

EXTINGUISHING MEDIA: Carbon dioxide, water, dry chemical

SPECIAL FIRE FIGHTING PROCEDURES: Avoid extinguishing methods which may generate dust clouds. Water stream can disperse dust in air, producing a fire hazard and possible explosion hazard if exposed to ignition source.

UNUSUAL FIRE & EXPLOSION HAZARDS: Polymer dust is combustible. The explosive limits of the polymer particles suspended in air are approximately those of coal dust. Firefighters should wear self-contained breathing apparatus.

SECTION V: HEALTH HAZARD DATA

INGESTION: Rinse mouth out with water. Ingestion of small quantities of this material under normal circumstances would not cause harmful effect.